

# Break Through Traditional Swab Limitations

The Rhinostics GrooveSwab™ Nasopharyngeal Swab ushers in a new paradigm in comfortable, convenient swab sample collection and processing workflows.

## Developed to Enhance Performance

Specialized hydrophobic biopolymer design includes a snap point to suit most tube configurations and enables full, rapid particle release during processing.

## Flexible for Your Needs

GrooveSwabs are individually packaged/sterilized and available in flocked or unflocked configurations and with or without a collection tube.



## Patient Comfort

The proprietary honey dipper swab head design is more comfortable than traditional swabs during sample collection.

## Dry Transport

No viral transport media (VTM) is used, which enables 30x increase in collection yield, reduced reagent costs and biohazard exposure, and less need for time-consuming concentration steps compared to traditional VTM-based swab transport methods.

## Ordering Information

- **GW-UF100001:** Rhinostics GrooveSwab Unflocked
- **GW-F100001:** Rhinostics GrooveSwab Flocked

## Upgrade Your Nasopharyngeal (NP) Swabs

It's time to replace your old NP swab. Learn more, order a sample, or request a quote at [www.rhinostics.com](http://www.rhinostics.com), or email:

✉ [sales@rhinostics.com](mailto:sales@rhinostics.com)