

AlteCap™ Switch: Handheld Screw Cap Decapper for RHINOstic® Automated Swab



Features:

- Semi-automatic, handheld device for decapping and capping 96 devices from SBS rack format sample tubes.
- Compatible with RHINOstic® Automated Swabs and Collection Devices and other screw cap tubes in SBS rack format.
- Can be used with racks of tubes in the portrait or landscape position.
- Easy to use, intuitive touch button controls.
- Light and portable.
- Internal cap driver control ensures each device is perfectly capped.
- Recap tubes with devices or caps. Apply fresh caps from a cap carrier.
- Multi-use stand is convenient for storage, USB charging, and holding caps for recapping.
- Cap driver cassettes are available for different types and brands of screw cap tubes in 24, 48, and 96 SBS rack formats.

Specifications:

Dimensions (LxWxH)	190 x 125 x 820 mm
	7.5 x 4.9 x 32 in
Weight	726 g alone/ 1.5 kg with stand
	1.6 lbs. alone/ 3.4 lbs with stand
Electrical	100-240V/50-60Hz
Operating Temperature Range	20-22° C/68-72° F
Connection	Micro-USB Charging Cable
Additional Supported Labware	SBS format labware

Compatible Rhinostics Products:

Part Number	Product Name
RH-S000001	RHINOstic Automated Swab Device, Sterile, Clear, Unmarked Bulk Packaged Capped Tubes
LV-2D1D002	RHINOstic Automated Swab Device, Sterile, Clear Uncapped 1D/2D/Visible Barcoded Tube
LV-2D1D003	RHINOstic Automated Swab Device, Sterile, Black Uncapped 1D/2D/Visible Barcoded Tube
VS-S000011	VERIstic Collection Device with Bulk Capped Tubes
LV-RK96001	96-Well Rack, Plastic
RR-100002	RHINOrac Aluminum Rack
RRT-A00001/2	RHINOrest/RHINOrest2 Aluminum Cap Holder, Assembled (block and base)

**Rhinostics Inc.**

411 Waverley Oaks Road, Suite 118

Waltham, MA 02452

(877) 746-6789

(877) RINOSTX

**RHINOSTICS****Rhinostics Inc.**

quality@rhinostics.com

www.rhinostics.com