

PRODUCT SPECIFICATIONS**Rhinostics GrooveSwab® Nasopharyngeal Swab****Description:**

Innovative and proprietary, fully molded flexible sterile nasopharyngeal swab with stacked-ring head for more comfortable nasopharyngeal collection, faster elution, and increased assay sensitivity. Increase testing rate exponentially over typical manual process flow. Patents pending.

Physical Properties:

Material: Polypropylene

Color: Clear Overall length: 6.0 inches (15.0 cm)

Diameter of collection zone: 0.12 inches (3 mm)

Length of collection zone: 0.12 inches (3 mm)

Handle diameter: 0.1 inches (2.5 mm)

Snap-Point: 3.0 inches (76.2 mm) from tip of the swab

Qualities:

Medical grade quality

Highly flexible comfortable tip

Break point for collection tip transport

High article collection capacity

Fits various collection tubes, including 76 x 13 mm tubes

Manufactured in China

Applications:

Nasopharyngeal collection by medical personnel

Contact Information:**Rhinostics Inc.**

411 Waverley Oaks Road, Suite 118

Waltham, MA 02452

(877) 746-6789

(877) RINOSTX

**Rhinostics Inc.**

quality@rhinostics.com

www.rhinostics.com