

## PRODUCT SPECIFICATIONS

# RHINOstic™ Nasal Swab Collection Device

### Description:

Innovative and proprietary, fully molded flexible sterile nasal swab with integral de-capper cap. Increase testing rate exponentially over typical manual process flow. Patents pending.

### Physical Properties:

Material: Polypropylene

Color: Clear

Overall length: 1.92 inches (48.7mm)

Diameter of collection zone: .12 inches (3mm)

Length of collection zone: .61 inches (15.4mm)

Cap diameter: .34 inches (8.6mm)

### Qualities:

Medical grade quality

Highly flexible comfortable tip

High article collection capacity

Fits standard 1 ml cryotube

Made in the USA

### Applications:

Anterior nares nasal swab collection

### Contact Information:



#### Rhinostics Inc.

411 Waverley Oaks Road, Suite 118

Waltham, MA 02452

(877) 746-6789

(877) RINOSTX



#### Rhinostics Inc.

[quality@rhinostics.com](mailto:quality@rhinostics.com)

[www.rhinostics.com](http://www.rhinostics.com)