

PRODUCT SPECIFICATIONS**Rhinostics Standard Nasal Swab****Description:**

Innovative and proprietary, fully molded flexible sterile nasal swab with break-off tip. Patents Pending.

Physical Properties:

Material: Polypropylene

Color: Blue

Overall length: 3.0 inches (76.0 mm)

Diameter of collection zone: Tapered 0.17 to 0.20 inches (4.3 to 5.0 mm)

Length of collection zone: 0.68 inches (17.27 mm)

Qualities:

Medical grade quality

Flexible comfortable tip

Break point for collection tip transport

High article collection capacity

Fits standard 1 ml cryotube

Made in the USA

Applications:

Anterior nares nasal collection

Contact Information:**Rhinostics Inc.**

411 Waverley Oaks Road, Suite 118

Waltham, MA 02452

(877) 746-6789

(877) RINOSTX

**RHINOSTICS**

Rhinostics Inc.

quality@rhinostics.com

www.rhinostics.com