

PRODUCT SPECIFICATIONS**VERIstic™ Collection Device****Description:**

Unique 3-sided capillary device with integral decapping cap. Combines wicking action to collect blood sample with automated processing to save time and reduce labor. Patent pending.

Physical Properties:

Material: Polypropylene

Color: Clear

Overall length: 1.9 inches (48 mm)

Cap diameter: 0.34 inches (8.6mm)

Qualities:

Medical grade quality

Three-sided capillary design

Fits standard 1 ml cryotube

Made in the USA

Applications:

Microsampling Blood Collection

Contact Information:**RhinoStics Inc.**

411 Waverley Oaks Road, Suite 118

Waltham, MA 02452

(877) 746-6789

(877) RINOSTX

**RhinoStics Inc.**

quality@rhinostics.com

www.rhinostics.com