

AlteCap™ Solo: Handheld Capper/Decapper for RHINOstic® Automated Swab



Features:

- Quickly cap, or recap, single tubes.
- Compatible with RHINOstic® Automated Swabs and Collection Devices and other screw cap tubes in SBS rack format.
- Easy to use, intuitive touch button controls.
- Light and portable.
- Internal cap driver control ensures each tube and automated device is perfectly capped.
- Recap tubes with devices or caps. Apply fresh caps from a cap carrier.
- Multi-use stand is convenient for storage, USB charging, and holding caps for recapping.

Specifications:

Dimensions (WxH)	25 x 190 mm
	.98 x 7.5 in
Weight	156 g alone/383 g with stand
	5.5 oz alone/13.5 oz with stand
Electrical	100-240V/50-60Hz
Operating Temperature Range	20-22°C/68-72°F
Connection	Micro-USB Charging Cable
Additional Supported Labware	SBS format labware

Compatible RhinoStics Products:

Part Number	Product Name
RH-S000001	RHINOstic Automated Swab Device, Sterile, Clear, Unmarked Bulk Packaged Capped Tubes
LV-2D1D002	RHINOstic Automated Swab Device, Sterile, Clear Uncapped 1D/2D/Visible Barcoded Tube
LV-2D1D003	RHINOstic Automated Swab Device, Sterile, Black Uncapped 1D/2D/Visible Barcoded Tube
VS-S000011	VERIstic Collection Device with Bulk Capped Tubes
LV-RK96001	96-Well Rack, Plastic
RR-100002	RHINOrac Aluminum Rack
RRT-A00001/2	RHINOrest/RHINOrest2, Aluminum Cap Holder, Assembled (block and base)



RhinoStics Inc.

411 Waverley Oaks Road, Suite 118

Waltham, MA 02452

(877) 746-6789

(877) RINOSTX



RhinoStics Inc.

quality@rhinostics.com

www.rhinostics.com